# Micro Prep

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
Treating with pinpoint precision and minimal invasive in the narrowest space is the declared goal of modern dentistry.

12 new Micro-Prep-Instruments are our contribution to this.

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
Specially designed instruments are required for a selective, minimal invasive preparation. A minimum working part diameter, application-oriented shapes, two diamond grit sizes, long thin necks and a perfect concentricity are the essential characteristics of such instruments.

The long, thin neck construction enables perfect optimal control of preparation deep inside the dentin, even in hard-to-reach places. Specific work can be done without opening the enamel to wide, which can be the case using the same working part diameters with usual neck construction.

Hardened special steel enables safe working at a maximum compressive force of 1-2 N that gives a new softer grinding feeling. The instruments with FG-shank come in six different shapes, each in medium and fine diamond grit for the purpose of defined abrasion and optimal surface quality.

We recommend a working speed of 160.000 RPM and spray cooling at least (50ml/min). All current disinfection, cleaning and sterilization methods are acceptable. CE-mark on every package.

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Interested in an offer? Just fill in this form with the required number of pieces and fax it to BUSCH!

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
Micro-Prep-Instruments also available as a set! (1 piece per type on STERI-SAFE stands)

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