



## NON STERILE SLEEVES WITH ADHESIVE

Polyethylene protection sleeves cover handpieces, tubes, turbines and contra-angles.

- Convenient: packaging makes them easy to tear-off
- Easy to set up
- Preserve the tactile sensation of the surgeon while handling equipment: manufactured in a flexible and adherent material

② Non sterile single use



Box of 250 sleeves  
Ref. HY-066

The sleeve's dimensions allow the protection of tubing 90cm in length and 8.5cm in width.

The sleeves contain an adhesive which provides good grip on the tubing.

## 45MM WIDE SLEEVES FOR TURBINES

Length: 45cm, 1.5cm wide conical end

② Non sterile single use



Roll of 250 sleeves  
Ref. HY-060



## 55MM WIDE SLEEVES FOR CONTRA-ANGLES, HANDPIECES AND WATER/AIR SYRINGES

Length: 45cm, 1.5cm wide conical end

② Non sterile single use



Roll of 250 sleeves  
Ref. HY-062

## STERILE SLEEVES

Sterile sleeves provide excellent protection for tubing while guaranteeing the best sanitary conditions.

This disposable PVC sleeve, specifically designed for implantologists, is both:

- **Easy to put on:** on account of its user-friendly packaging. An inner card enables the operator to extend the sleeve without having to touch the tube
- **Easy to secure** using the latex-free elastic bands supplied

Dimensions correspond to standard tubing size: 120cm long x 8cm wide.



② Sterile single use



1 single sleeve  
Ref. HY-138



Sets of 2 sleeves  
Ref. HY-139